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J. K. 73 Jahre alt

Evidenz

Initial Experience and Results With a New Adjustable Transobturator Male System for the Treatment of Stress Urinary Incontinence

Purpose: We report on our initial experience in terms of efficacy and safety with a new, self-anchoring adjustable transobturator male system (A.M.I.® ATOMS System) for the treatment of male stress urinary incontinence after prostate surgery.

Materials and Methods: In this prospective, nonrandomized single center study conducted between March and December 2009, patients with stress urinary incontinence secondary to prostatic surgery were treated with the ATOMS device. Urethroscopy, filling and voiding cystometry were performed pre-operatively for all patients. In addition, incontinence symptoms were assessed, and a physical examination, 24-hour pad test and 24-hour pad count were performed before and after surgery. **Results:** A total of 38 patients were included in the study (36 after radical prostatectomy, 2 after benign prostatic hyperplasia surgery). No intraoperative complications occurred. Mean number of adjustments during followup was 3.97 (range 0 to 9). At a mean followup of 16.9 months (range 13 to 21) the overall success rate was 84.2%. Of the successful cases 60.5% were considered dry (0 to 1 pad and less than 15 ml/24-hour pad test) and 23.7% improved (more than 1 pad per 24 hours but more than 50% decrease in pad use and less than 100 ml per 24-hour pad test). In 15.8% of the patients the treatment was considered to have failed (more than 2 pads daily and greater than 100 ml on 24-hour pad test). **Conclusions:** The treatment of male stress urinary incontinence with the ATOMS is safe and effective. It is an excellent first or second line treatment for mild to moderate male stress urinary incontinence, even after external irradiation. The option of longterm, minimally invasive adjustment to respond to patient needs is a significant advantage of this new implant.

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M.R. Hoda - Der Urologe A July 2012; pp. 1-9

Erfolgsrate 93,8%

Early results of a European multicentre experience with a new self-anchoring adjustable transobturator system for treatment of stress urinary incontinence in men

Objective: To report our experience with a new self-anchoring adjustable transobturator male system (ATOMS®; AMI, Austria) for the treatment of stress urinary incontinence (SUI) in men. **Patients and Methods:** A total of 99 men, in a number of centres, were treated for SUI with the new ATOMS® device. The device was implanted in all patients using an outside-in technique by passing the obturator foramen and anchoring the device to the inferior pubic ramus. The titanium port was placed s.c. on the left symphysis region. Adjustments were performed via port access. Postoperative evaluation consisted of physical examination, 24-h pad test, and 24 h-pad count. Preoperatively and at 6-month follow-up, patients completed a validated quality-of-life questionnaire. Two-way ANOVA was used to analyse changes over time. Within-group effects for time were tested using post hoc Dunnett's contrasts of baseline values vs subsequent measurements. **Results:** The most common indication was SUI after radical prostatectomy (92.9%). Failure of previous surgeries was present in 34.3% patients and 31.3% patients had undergone secondary radiation. The mean (SD; range) surgery time was 47 (13.8; 29–112) min. Temporary urinary retention occurred in two patients (2%) and transient perineal/scrotal dysaesthesia or pain was reported by 68 patients (68.7%) and resolved after 3–4 weeks of non-opioid analgesics. There were four (4%) cases of wound infection at the site of the titanium port leading to explantation. No urethral or bladder injuries related to the device or erosions occurred. The mean (SD; range) number of adjustments to reach the desired result (dryness, improvement and/or patient satisfaction) was 3.8 (1.3; 1–6). After a mean (SD; range) follow-up time of 17.8 (1.6; 12–33) months, the overall success rate was 92% and the mean pad use decreased from 7.1 to 1.3 pads/24 h (P < 0.001). Overall, 63% were considered dry and 29% were improved. **Conclusion:** Treatment of male SUI with this self-anchored adjustable system is safe and effective.

M.R. Hoda / G. Primus - BJU Int 2013 Feb;111(2):296-303

Patientenzufriedenheit 92%

Five-year experience with the adjustable transobturator male system for the treatment of male stress urinary incontinence: a single-center evaluation

Background: We report on our 5-year experience with the adjustable transobturator male system (ATOMS®, A.M.I., Feldkirch, Austria).

Methods: Between 10-2009 and 10-2014, 54 patients received an ATOMS. The mean follow-up of this retrospective observational trial was 27.5 ± 18.4 (2.3–59) months. Within each follow-up, the following were evaluated: micturition protocol, 24-h pad count, uroflowmetry and residual volume. Statistical analysis was performed with Sigma-Plot® 11.0, p < 0.05 considered as significant. **Results:** Stress urinary incontinence (SUI) I°, II° and III° was seen in 1 (1.9 %), 16 (29.6 %) and 37 patients (68.5 %), respectively. In summary, 48.1 % of the patients became "dry" (0-"safty pad"/day), while 29.6 % achieved at least an "improvement" of about more than 50 % (1–2 pads/day), which corresponds to an overall success rate of 77.7 %. The mean number of pads/day decreased from 7.7 to 1.6. Regarding the initial degree of SUI, patients with mild or moderate incontinence had a significantly better outcome (p = 0.002, 95 % CI 0.9066 to 2.760). Postoperative complications were scaled according to the Clavien classification, in which we have seen 4 grade I-, 1 grade IIIa- and 9 grade IIIb-complications (overall 25.9 %). The evaluation of quality of life by ICIQ-SF showed a significant improvement (p = 0.0001, 95 % CI -14.56 to -11.75). **Conclusion:** The treatment of male SUI using the ATOMS incontinence system achieved the best results in patients with mild and moderate incontinence. For severe incontinent patients, the system represents an efficient alternative.

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International patent filed / pending / granted

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A.M.I. GmbH
A.M.I. Deutschland GmbH
FUMEDICA AG