

StatLock® CV Plus Stabilization Device

“Suture-Free” Stabilization for Central Venous Catheters

Minimize Overall Central Venous Catheter Complications.^{1,4,5}

Maximize Central Venous Catheter Patient Comfort.

- Designed to reduce the potential complications associated with suture securement^{1,4,5}
- Designed to improve patient comfort and safety by eliminating the need to secure catheters with suture
- Prevents needlestick injuries associated with catheter securement^{2,4,5}
- Provides solutions that follow OSHA's recommendations for sutureless securement²

***An Evidence-Based
Stabilization Device
That Meets The
2006 INS Standards³***



StatLock®

Safety. Security. Savings.

StatLock® CV Plus Stabilization Device

Clinical studies show that StatLock® devices

- Eliminate accidental needlesticks associated with catheter securement.^{2,4,5}
- Reduce catheter complications.^{1,4,5}

StatLock® Device Retainer

- Secures catheter in place.
- Accommodates a wide variety of central line catheters.

StatLock® Device Anchor Pad

- Faster and easier to apply than suture.^{4,5}
- Breathable latex-free tricot fabric provides added catheter security and patient comfort.
- Tincture of Benzoin included to enhance securement.

StatLock® Device Pigtail Retainer

- Ideal for organizing and separating catheter extensions.
- Reduces tension on catheter extensions.



Ordering Information

Reference Number	Description	QTY/Box
CV0220	StatLock® CV Plus Device, Crescent Tricot, Sliding Post, Pigtail Retainer	25
CV0223	StatLock® CV Plus Device, Crescent Tricot, Fixed Post, Pigtail Retainer	25
CV0224	StatLock® CV Plus Device, Butterfly Tricot, Sliding Post, Pigtail Retainer	25
CV0226	StatLock® CV Plus Device, Crescent Tricot, Sliding Post, No Pigtail Retainer	25
CV0321	StatLock® CV Ultra Device, Tricot, Fixed Post, Tethered Locking Caps, No Pigtail Retainer	25
VCDPP	StatLock® CV Plus Pediatric Device, Foam, Sliding Post, Pigtail Retainer	25
VCDT	StatLock® CV Plus Pediatric Device, Foam, Fixed Post, No Pigtail Retainer	50
VCDPS	StatLock® Pigtail Retainer Device	50

StatLock® devices are sterile and latex-free.



BARD MEDICAL DIVISION

C.R. Bard, Inc.
8195 Industrial Blvd.
Covington, GA 30014

Customer Service 800.526.4455
Fax 800.852.1339

www.statlock.com

Please consult product inserts and labels for any indications, contraindications, hazards, warnings, cautions and directions for use.

Bard and StatLock are registered trademarks of C.R. Bard, Inc. or an affiliate.

Venetec International, Inc. is a subsidiary of C.R. Bard, Inc.

©Copyright 2006 C.R. Bard, Inc. All rights reserved.

1. Wood D, Bowe-Geddes LA. A comparative retrospective analysis of two securement techniques for peripherally inserted central catheters (PICC) and Midlines in the homecare setting. JVAD 1997 Vol 2, No 3:11-16.
2. OSHA Fact Sheet #7. Securing Medical Catheters. 2004.
3. INS Standards of Practice, Section 43B, 2006.
4. Yamamoto AJ, et. al. Sutureless Securement Device Reduces Complications of Peripherally Inserted Central Venous Catheters. J Vasc Interv Radiol 2002 13:77-81
5. Schears G, et al. Spotlight on Solutions, JCAHO/NPSF 2000.