

L-HOOK
electrode

VECTEC®
SINGLE USE
LAPAROSCOPIC
DEVICES





Ref. : CAB 3000
Disposable monopolar cable

disposable laparoscopic monopolar L-hook dissector

Reference	Working length	Diameter	Active part
H5 330	330 mm	5 mm	axial: 5 mm radial: 5 mm
H5 330L	330 mm	5 mm	radial: 5 mm

The L-hook electrode with blunt tip ensures fine dissection.

Complete and even sheathing of the instrument guarantees optimum electrical insulation and comfortable grip during the surgical procedure.

The integrated hook connector fits any monopolar cable plug which conforms to current safety standards.

Packaging : box of 25 individually packaged units

MD Class IIb
Sterilized by ETO
See instructions for use



REF : H5 330



REF : H5 330L

Manufactured by:

Péters Surgical
42 rue Benoit Frachon
93013 Bobigny Cedex - FRANCE

Tel: +33 (0)1 48 10 62 62
Fax: +33 (0)1 48 91 22 99
commercial@peters-surgical.com
www.peters-surgical.com