



Made in  
the **EU**

99.  
9999%

EFFECTIVE **UV-C**  
DISINFECTION  
TECHNOLOGY

**steri**pro



0.0001%

tolerance to  
any kind of infection

steripro



“My mother went to hospital to have kidney stones removed. The operation was a success, but she developed sepsis unexpectedly, which ended fatally.”

**Jane, 26**



“My grandfather had to have a routine blood vessel operation. During his recovery in the hospital, he got COVID-19 and died three weeks later as a result of the infection. I wish hospitals had better prevention procedures.”

**Marco, 31**



“My sister underwent brain surgery, which went well. But during her recovery in hospital she got an infection that led to her losing her sight.”

**Anna, 54**

# Do you want safer hospitals for yourself and your loved ones?





Hospital grade  
**UV-C**  
disinfection  
device



**1 EUR /**  
disinfection



# Let's prevent hospital associated infections

- Manual cleaning and chemical disinfection leave less than 50% of hospital room surfaces adequately cleaned when chemical germicides are used\*
- In the EU, 110,000 people die each year from infections, which is 4 times more than the number of road traffic fatalities\*\*
- Each year in the EU, over 4 million patients acquire a healthcare-associated infection (E. coli, MRSA, VRE, COVID-19)\*\*\*
- Because viruses and bacteria are invisible to the human eye, they remain on surfaces even after manual cleaning\*\*\*


## Sources:

\* D. J. Weber and W. A. Rutala. Understanding and Preventing Transmission of Healthcare-Associated Pathogens Due to the Contaminated Hospital Environment. Infection control and Hospital Epidemiology. May 2013, vol. 34, no. 5.

\*\* 2019 road safety statistics: what is behind the figures?  
([https://ec.europa.eu/commission/presscorner/detail/en/qanda\\_20\\_1004](https://ec.europa.eu/commission/presscorner/detail/en/qanda_20_1004))

\*\*\* Health at a Glance: Europe 2016  
([https://www.oecd-ilibrary.org/social-issues-migration-health/health-at-a-glance-europe-2016/healthcare-associated-infections\\_health\\_glance\\_eur-2016-46-en;jsessionid=9C4QCdHRW6vFoWKLDOC672rr.jp-10-240-5-167](https://www.oecd-ilibrary.org/social-issues-migration-health/health-at-a-glance-europe-2016/healthcare-associated-infections_health_glance_eur-2016-46-en;jsessionid=9C4QCdHRW6vFoWKLDOC672rr.jp-10-240-5-167))

# Monitor. Evaluate. Disinfect.




SteriPro is a smart device that offers a smart solution to your disinfection challenges. It operates as a seamless ecosystem, through the SPOT portal, which gives you total control of device utilisation and status tracking. **SteriPro enables up to 70 quick disinfections per day (5 to 15 min cycles) with validated 99.9999% efficiency and therefore enables hospitals to save a lot of time.**




# Certified and validated disinfection


The validated and certified SteriPro's technology is proven to be exceptionally efficient at disinfection, thanks to its ultraviolet light. Mercury vapor lamps radiate about 95% of their energy at a wavelength 253.7 nm, which is close to the DNA absorption peak. This means SteriPro has a high germicidal effectiveness against microbes, rendering the organisms sterile in no time. **Confirmed and tested by independent laboratory.**



**MRSA**



**SARS-COV-2**



**VRE**

# Disinfect more (smartly)



Efficiency 99.9999%  
(log 6 reduction)



Certificated  
and validated



Disinfection cycles  
5 to 15min for 63 m<sup>2</sup>



Upgraded manual cleaning  
processes with contactless  
disinfection



Reaches all surfaces -  
even shadow areas in  
a radius of 12m



The device can work  
24/7, no charging, up to  
70 disinfections per day



Developed in  
partnership with  
hospitals



Highest level  
of safety



€1 = 1 disinfection



Durability guaranteed.  
Made in the EU



**steripro**

**Pay more attention to hygiene and  
protect yourself and your patients.**

Join hospitals from more than 30 countries around the world.



[www.uvc-solutions.com](http://www.uvc-solutions.com)



[info@uvc-solutions.com](mailto:info@uvc-solutions.com)



steripro

POWERED BY  
**PHILIPS UV-C LAMPS**

