

A.M.I. Product Catalog International

Issue 03/2025



Medical progress
made by A.M.I.

**Empowering Innovation.
Empowering Lives.**



The A.M.I. vision is clearly defined: With innovative products and surgical techniques, we improve medical treatment for people around the world. In the surgical disciplines of coloproctology, urology and urogynecology, we aim to be one of the leading medical technology companies and partners.

We inspire our customers with competence and trendsetting innovations. Inspiration for new ideas comes from both the market and modern findings from research and science. In this way, we always find the balance between decisive factors such as cost and time pressure, ergonomics, user-friendliness and, of course, above all, the well-being of patients. Through perfection, we generate the highest quality and patient safety day after day.

Holding ISO 13485 certification, we comply with globally recognized standards for quality management in medical devices.

The products contained in this overview meet the requirements of the Medical Device Directive 93/42 ECC or the Regulation (EU) 2017/745. The certificate covers all activities – from design and manufacturing to sales and shipment of products. For further information or queries, please contact us:

A.M.I. Agency for Medical Innovations GmbH
Im Letten 1
A-6800 Feldkirch
T: +43 (0) 5522 90505-0
E: customerservice@ami.at

We are looking forward to your request!

The A.M.I. Team






- Hemorrhoids
- Fecal Incontinence
- Anal Abscesses & Fistulas

Coloproctology



TRILOGY (Wi-3 HAL-RAR Unit)

Order Code		Product description	Included in delivery
TRI2010		TRILOGY Unit (Wi-3 HAL-RAR Unit) The Wi-3 HAL-RAR System (HAL - Hemorrhoidal Artery Ligation, RAR - Recto Anal Repair) is a Bluetooth-enabled, pulsed Doppler, ultrasound system intended to be used for detecting hemorrhoidal arteries, to manually ligate these arteries (HAL), and subsequently carry out a mucopexy of hemorrhoidal tissue prolapse (RAR) if necessary	Delivered non-sterile, steam autoclavable 1 system consists of: <ul style="list-style-type: none"> • Doppler-ultrasound electronic unit with integrated bluetooth submitter • TRI2020 Wi-3 Battery Cap • TRI2030 Wi-3 Fixation nut • TRI2040 Wi-3 Battery Set Multi-use (50 applications) Batteries: 2 units, AA NiMH, 2100 mAh
TRI2110		Speaker (Wi-3 Speaker) A non-medical, non-sterile and Bluetooth-enabled accessory component of the TRILOGY Unit (Wi-3 HAL-RAR Unit)	1 unit, supplied with USB cable
TRI2210		Adapter A non-medical, non-sterile USB adaptor for connecting the speaker Manufacturer Frimo FOX12-X11-USB	1 unit, delivered non-sterile
TRI2220		Connector Europe (Type C) A non-medical, non-sterile, replaceable connector for connection to adaptor, Manufacturer Frimo	1 unit, delivered non-sterile
TRI2020		Wi-3 Battery Cap Reusable battery cap for the TRILOGY Unit (Wi-3 HAL-RAR Unit)	1 unit, delivered non-sterile, steam autoclavable
TRI2030		Wi-3 Fixation Nut A reusable connecting element which enables a form-fit connection between TRILOGY Unit (Wi-3 HAL-RAR Unit) and probe	1 unit, delivered non-sterile, steam autoclavable


TRILOGY (Wi-3 HAL-RAR Unit)

Order Code		Product description	Included in delivery
TRI2070		Wi-3 Arm Reusable accessory to comfortably hold the Wi-3 HAL-RAR system in position	1 unit, delivered non-sterile, steam autoclavable
TRI2040		Wi-3 Battery Set 1 package consisting of 2 rechargeable batteries (AA, NiMH, 2100 mAh)	2 pieces, delivered non-sterile
102372		Wi-3 Battery Charger Suitable for recharging up to 4 batteries	1 unit, delivered non-sterile
AHN 006		A.M.I. HAL Needleholder Reusable needle holder for use during hemorrhoidal artery ligation For use with Wi-3 HAL-RAR system (TRI2010)	1 instrument, delivered non-sterile, steam autoclavable
AHK 007		A.M.I. HAL Knotpusher Reusable knot pusher for placing knots during hemorrhoidal artery ligation For use with Wi-3 HAL-RAR system (TRI2010)	1 instrument, delivered non-sterile, steam autoclavable
AHAL70		A.M.I. HAL Suture Suture material for HAL and RAR procedures.	1 box, delivered sterile 36 x sutures, 5/8 needle absorbable, 2/0, 75cm

+ TRILOGY (Wi-3 HAL-RAR Unit) Probes

Order Code		Product description	Included in delivery
RAR2081		RAR Flexi Probe Single-use probe with sleeve for a variable opening window during detection, ligation and mucopexy of hemorrhoidal arteries	1 box, delivered sterile 5 x RAR Flexi probe
RAR2181		HAL-RAR Procedure Kit – Flexi	1 kit, delivered sterile 1 x RAR Flexi probe 8 x A.M.I. HAL Suture 1 x Single-use knotpusher designed especially to fit the A.M.I. probes

Soft Anal Band

Order Code		Product	Included in delivery
AAB 300		Soft Anal Band Long-term subcutaneous implant to control the sphincter function for the treatment of fecal incontinence in adults	1 system, delivered sterile <ul style="list-style-type: none"> • A.M.I. Soft Anal Band Ring (cuff) in two different sizes M-L • A.M.I. Soft Anal Band Valve • A.M.I. Soft Anal Band Activator M-L • A.M.I. Safe Puncture Port Low Profile • A.M.I. Port Needle • Ring Connector • Port Connector • Patient card, Patient lot labels





Comfort Drain

Order Code	Product description	Included in delivery
CDS2001	Comfort Drain A flexible, silicone, elastic drain with a knotless closure mechanism to facilitate drainage of perianal fistulas	1 box, delivered sterile 10 x set
CDS2010	Comfort Drain Grasper Reusable instrument to facilitate the closure of the Comfort Drain CDS2001	1 box, non-sterile, steam autoclavable 1 x Comfort Drain Grasper

FiXcision

Order Code	Product	
FIX2001	FiXcision A surgical invasive multi-part cutting instrument for anal fistula treatment. FiXcision consists of the following components: <ul style="list-style-type: none">• Probe• Guide• Base Plate• Circular Cutter	1 box, delivered sterile 10 x set

Coloproctology Demo Products

Order Code		Product description	Included in delivery
TRI2010-DEMO		TRILOGY Demo Wi-3 HAL-RAR Unit DEMO for unlimited use	1 unit, delivered non-sterile
RAR2081-DEMO		RAR Flexi Probe DEMO	1 set (probe and sleeve) delivered non-sterile
TRI-HOLDER		TRILOGY Demo Holder Plexiglass holder for presentation of TRILOGY. Can also be used as a training model with clamped foam. (TRILOGY not included)	1 holder, delivered non-sterile
FIS-DEMO		Fistula Demo Case Set of 5 Comfort Drains 1 Comfort Drain Grasper 1 FiXcision 1 forceps	1 demo case All products are non-sterile

All demo products are strictly for demonstration purposes only and not intended for use in the operating theatre!

Treatment of Pelvic Organ Prolapse (POP)

Single-incision vaginal mesh procedure



- InGYNious
- BSC Mesh

Laparoscopic mesh procedure

- EndoGYNious
- PelviGYNious
- ProGYNious

Instruments

- i-Stitch



Treatment of Stress Urinary Incontinence (SUI)

Retropubic and trans-obturator mid-urethral sling procedure

- sensiTVT
- sensiTVT-A
- A.M.I. TVA Sling
- A.M.I. TOA Sling
- A.M.I. Multi Purpose Sling
- A.M.I. Tunneller

Contained laparoscopic power morcellation





- i-Cut
- More-Cell-Safe



Solutions for
Urogynecology

Treatment of Pelvic Organ Prolapse (POP)





Single-incision vaginal mesh procedure – InGYNious

Order Code		Product description	Included in delivery
IGY5961		InGYNious D A S * Ultralight, polypropylene mesh with i-Stitch loading units for 3-level support of anterior vaginal wall prolapse	1 box, all items delivered sterile: 1 x InGYNious D A S (Standard) 6 x i-Stitch loading unit PET 0 W (from IST1031) • non-absorbable, polyester suture, USP 0
IGY5951		InGYNious D A L * Ultralight, polypropylene mesh with i-Stitch loading units for 3-level support of anterior vaginal wall prolapse	1 box, all items delivered sterile: 1 x InGYNious D A L (Large) 6 x i-Stitch loading unit PET 0 W (from IST1031) • non-absorbable, polyester suture, USP 0
IGY5981		InGYNious D P S * Ultralight, polypropylene mesh with i-Stitch loading units for 3-level support of posterior vaginal wall prolapse	1 box, all items delivered sterile: 1 x InGYNious D P S (Standard) 4 x i-Stitch loading unit PET 0 W (from IST1031) • non-absorbable, polyester suture, USP 0
IGY5971		InGYNious D P L * Ultralight, polypropylene mesh with i-Stitch loading units for 3-level support of posterior vaginal wall prolapse	1 box, all items delivered sterile: 1 x InGYNious D P L (Large) 4 x i-Stitch loading unit PET 0 W (from IST1031) • non-absorbable, polyester suture, USP 0

* D = Direct

A = Anterior
S = StandardP = Posterior
L = Large

PP = Polypropylen



Order Code		Product description	Included in delivery
IGY5561		InGYNious D A S-PP * Ultralight, polypropylene mesh with i-Stitch loading units for 3-level support of anterior vaginal wall prolapse	1 box, all items delivered sterile: 1 x InGYNious D A S-PP (Standard) 6 x i-Stitch loading unit PP 0 (from IST1011) • non-absorbable, polypropylene suture, USP 0
IGY5551		InGYNious D A L-PP * Ultralight, polypropylene mesh with i-Stitch loading units for 3-level support of anterior vaginal wall prolapse	1 box, all items delivered sterile: 1 x InGYNious D A L-PP (Large) 6 x i-Stitch loading unit PP 0 (from IST1011) • non-absorbable, polypropylene suture, USP 0
IGY5571		InGYNious D P L-PP * Ultralight, polypropylene mesh with i-Stitch loading units for 3-level support of posterior vaginal wall prolapse	1 box, all items delivered sterile: 1 x InGYNious D P L-PP (Large) 4 x i-Stitch loading unit PP 0 (from IST1011) • non-absorbable, polypropylene suture, USP 0
IGY5581		InGYNious D P S-PP * Ultralight, polypropylene mesh with i-Stitch loading units for 3-level support of posterior vaginal wall prolapse	1 box, all items delivered sterile: 1 x InGYNious D P S-PP (Standard) 4 x i-Stitch loading unit PP 0 (from IST1011) • non-absorbable, polypropylene suture, USP 0

* D = Direct


A = Anterior
S = StandardP = Posterior
L = Large

PP = Polypropylene






Single-incision vaginal mesh procedure – InGYNious

Order Code		Product description	Included in delivery
IGY5921		InGYNious V Ultralight, polypropylene mesh with i-Stitch loading units for 3-level repair of vaginal vault prolapse	1 box, all items delivered sterile: 1 x InGYNious V 6 x i-Stitch Loading Unit PET 0 W (from IST1031) • non-absorbable, polyester suture, USP 0
IGY-BAR		InGYNious Bar Help tool for procedures for prolapse surgery with InGYNious	1 box, delivered non-sterile: 1 x InGYNious Bar




Single-incision vaginal mesh procedure – BSC Mesh

Order Code		Product description	Included in delivery
BSC5001		BSC Mesh Ultralight, polypropylene mesh designed to induce the formation of neo-ligaments by establishing symmetrical, bilateral suspension of the vaginal vault from the sacrospinous ligament	1 box, all items delivered sterile: 1 x BSC Mesh 2 x i-Stitch loading unit PP 0 (from IST1011) • non-absorbable, polypropylene suture, USP 0




i-Stitch

Order Code		Product description	Included in delivery
IST1010		i-Stitch Reusable instrument for pelvic floor reconstruction; for native tissue repair or synthetic mesh placement	1 box, delivered non-sterile, steam autoclavable: 1 x i-Stitch 1 x i-Stitch Cleaning Brush
IST1040		i-Stitch up Reusable instrument for pelvic floor reconstruction; for native tissue repair or synthetic mesh placement	1 box, delivered non-sterile, steam autoclavable: 1 x i-Stitch up 1 x i-Stitch Cleaning Brush
IST1011		i-Stitch Loading Unit PP 0 Disposable device consisting of a needle slider (incl. suture clamp), a cannulated needle, a suture dispenser, and a suture with a specially formed blunt, ball-shaped suture tip for the i-Stitch/i-Stitch up	1 box, delivered sterile: 5 x i-Stitch loading unit PP 0 • non-absorbable, monofil polypropylene suture, suture size USP 0 • Color: blue
IST1051		i-Stitch Loading Unit PDO 0 Disposable device consisting of a needle slider (incl. suture clamp), a cannulated needle, a suture dispenser, and a suture with a specially formed blunt, ball-shaped suture tip for the i-Stitch/i-Stitch up	1 box, delivered sterile: 5 x i-Stitch loading unit PDO 0 • absorbable, monofil polydioxanone PDO, suture size USP 0 • Color: blue
IST1021		i-Stitch Loading Unit PDO 2-0 Disposable device consisting of a needle slider (incl. suture clamp), a cannulated needle, a suture dispenser, and a suture with a specially formed blunt, ball-shaped suture tip for the i-Stitch/i-Stitch up	1 box, delivered sterile: 5 x i-Stitch loading unit PDO 2-0 • absorbable, monofil polydioxanone PDO, suture size USP 2-0 • Color: blue

i-Stitch






Order Code		Product description	Included in delivery
IST1031		i-Stitch Loading Unit PET 0 W Disposable device consisting of a needle slider (incl. suture clamp), a cannulated needle, a suture dispenser, and a suture with a specially formed blunt, ball-shaped suture tip for the i-Stitch/i-Stitch up	1 box, delivered sterile: 5 x i-Stitch loading unit PET 0 W <ul style="list-style-type: none"> • non-absorbable, braided polyester suture, suture size USP 0 • Color: white
IST1050		i-Stitch Rinsing Needle Miele Connector Rinsing needle for the i-Stitch instruments	1 box, delivered non-sterile: 1 x i-Stitch Rinsing Needle Miele Connector
IST1060		i-Stitch Rinsing Needle Luer-Lock Connector	1 box, delivered non-sterile: 1 x i-Stitch Rinsing Needle Luer-Lock Connector
IST1020		i-Stitch Cleaning Brush Brush for cleaning the i-Stitch instruments	1 box, delivered non-sterile: 2 x i-Stitch Cleaning Brushes

Laparoscopic mesh procedure: ProGYNious, EndoGYNious, PelviGYNious



Order Code		Product	Included in delivery
PFR5731		ProGYNious Pre-shaped mesh prosthesis for surgical treatment (laparotomy or laparoscopic approach) of pelvic floor prolapse	1 box, delivered sterile: 1 x ProGYNious (anterior + posterior mesh combined)
PFR5641		EndoGYNious Y-shaped mesh prosthesis for surgical treatment (laparotomy or laparoscopic approach) of pelvic floor prolapse	1 box, delivered sterile: 1 x EndoGYNious (anterior + posterior mesh combined)
PFR5651		PelviGYNious Pre-shaped mesh prosthesis for surgical treatment (laparotomy or laparoscopic approach) of pelvic floor prolapse	1 box, delivered sterile: 1 x PelviGYNious (anterior + posterior mesh, 2 pieces)

Treatment of Stress Urinary Incontinence (SUI)

Retropubic mid-urethral sling procedure Transobturator mid-urethral sling procedure



Order Code		Product description	Included in delivery
SUI5011		sensiTVT Sling of monofilament macropore polypropylene for treatment of female SUI with passively articulating joints and adaptable suburethral graft	1 box, delivered sterile: 1 x sensiTVT
SUI5021		sensiTVT-A Sling of monofilament macropore polypropylene for treatment of female SUI. Equipped with: • passively articulating joints and adaptable suburethral graft • sutures for post-operative tension adjustment	1 box, delivered sterile: 1 x sensiTVT-A
TVA5021		A.M.I. TVA Sling Adjustable sling for treatment of female stress urinary incontinence	1 box, delivered sterile: 1 x A.M.I. TVA Sling
TOA5121		A.M.I. TOA Sling Adjustable sling for treatment of female stress urinary incontinence	1 box, delivered sterile: 1 x A.M.I. TOA Sling
PFR5021		A.M.I. Multi Purpose Sling Sling for treatment of female stress urinary incontinence	1 box, delivered sterile: 1 x A.M.I. Multi Purpose Sling

Tunneller

Order Code		Product description	Included in delivery
TVA5030		A.M.I. TVA Tunneller Multiple use devices made of stainless steel for the placement of A.M.I. urogynecological sling implants (Retropubic technique)	1 box, delivered non-sterile, steam autoclavable: 1 x A.M.I. TVA Tunneller
TOA5130		A.M.I. TOA Tunneller Multiple use devices made of stainless steel for the placement of A.M.I. urogynecological sling implants (transobturator technique / outside-in / inside-out technique)	1 box, delivered non-sterile, steam autoclavable: 1 x Pair of A.M.I. TOA Tunneller

Contained Laparoscopic Power Morcellation

i-Cut

Order Code		Product description	Included in delivery
ICT5011		i-Cut Single-use laparoscopic morcellation device	1 box, delivered sterile: 1 x i-Cut incl. obturator
ICT5020		i-Cut Power Supply Typ C – Euro Reusable power supply unit for i-Cut ICT5011 Additionally available: ICT5030 i-Cut Power Supply Type G ICT5040 i-Cut Power Supply Type I ICT5050 i-Cut Power Supply Type A ICT5060 i-Cut Power Supply Type C - Korea	1 box, delivered non-sterile: 1 x i-Cut Power Supply (ICT5010) 1 x Country specific adapter (power cord)

More-Cell-Safe

Order Code		Product description	Included in delivery
MCS5111		More-Cell-Safe 0° System indicated for laparoscopic endoscopic tissue morcellation in gynecology. Including Visi-Shield for 0° optics to prevent camera contamination	1 box, delivered sterile: 10 x More-Cell-Bag 10 x Visi-Shield 0°
MCS5151		More-Cell-Safe 30° System indicated for laparoscopic endoscopic tissue morcellation in gynecology. Including Visi-Shield for 30° optics to prevent camera contamination	1 box, delivered sterile: 10 x More-Cell-Bag 10 x Visi-Shield 30°
MCS5131		More-Cell Bag System indicated for laparoscopic endoscopic tissue morcellation in gynecology. Made of elastic, cell-tight material, which can be insufflated with CO2	1 box, delivered sterile: 10 x More-Cell-Bag

Urogynaecology Demo Products


Order Code		Product description	Included in delivery
BSC-DEMO		BSC Presentation Plexiglass holder for presentation of BSC Mesh	1 presentation set, delivered non-sterile 1 x BSC mesh, non-sterile 2 x acrylic glasses with magnets 1 x metal base
PELVIGYNIUS-DEMO		PelviGYNious Presentation Plexiglass holder for presentation of PelviGYNious	1 presentation set, delivered non-sterile 1 x PelviGYNious mesh, non-sterile 2 x acrylic glasses with magnets 1 x metal base
PROGYNIUS-DEMO		ProGYNious Presentation Plexiglass holder for presentation of ProGYNious	1 presentation set, delivered non-sterile 1 x ProGYNious mesh, non-sterile 2 x acrylic glasses with magnets 1 x metal base
IGY-DEMO		InGYNious Presentation Plexiglass holder for presentation of InGYNious	1 presentation set, delivered non-sterile 1 x InGYNious Anterior mesh non-sterile 1 x InGYNious Posterior mesh non-sterile 2 x acrylic glasses with magnets 1 x metal base
IGY-BONE		InGYNious Pelvis Demonstration Model Pelvic bone with implanted InGYNious mesh	1 model, delivered non-sterile

All demo products are strictly for demonstration purposes only and not intended for use in the operating theatre!

Urogynaecology Demo Products

Order Code		Product description	Included in delivery
IST-DEMO		i Stitch Demo Foam Metal plate with foam to trigger the i-Stitch suture. Size: 15 x 15 x 2 cm	1 unit, delivered non-sterile
SUI-DEMO		sensiTVT Presentation Plexiglass holder with fixed sensiTVT. Shows the effect of sensiTVT design paraurethral	1 model, delivered non-sterile
MCS-DEMO		MCS Demo Bag MCS-Bag only no Visi-Shield included	1 unit, delivered non-sterile
MCS-HOLDER		MCS Demo Holder Metal ring for demonstration of MCS only No Product included	1 model, delivered non-sterile
ICT-DEMO		i-Cut Demo i-Cut Instrument only No Power Supply	1 unit, delivered non-sterile
ICT-HOLDER		i-Cut Demo Holder Metal plate for presentation of i-Cut	1 holder, delivered non-sterile (i-Cut not included)
110516		UROGYN Portfolio EN UroGyn folder includes all A.M.I. meshes and slings implants for taking out and showing to customer (no instruments)	1 portfolio, delivered non-sterile

All demo products are strictly for demonstration purposes only and not intended for use in the operating theatre!






Solutions for
Urology

ATOMS

Order Code		Product description	Included in delivery
ATS5041		ATOMS Adjustable system for treatment of male stress urinary incontinence System contains: <ul style="list-style-type: none"> • implant with soft, adjustable suburethral cushion • integrated, silicone-coated titanium port • non-coring port needles 	Single-use system, delivered sterile For the ATOMS surgery with ATS5041 the following product is needed: 1x TOA5130 TOA Tunneller or 1x ATS5010 ATOMS Tunneller
TOA5130		A.M.I. TOA Tunneller Reusable instrument of stainless steel for outside-in tunnelling. Treatment of male stress urinary incontinence with the A.M.I. ATOMS System	1 pair, helical (left and right), delivered non-sterile steam autoclavable
ATS5010		ATOMS Tunneller Reusable instrument of stainless steel for outside-in tunnelling. Treatment of male stress urinary incontinence with the A.M.I. ATOMS System	1 pair, helical (left and right), delivered non-sterile steam autoclavable
ATS5051		Scrotal Port for ATOMS Scrotal port with catheter, tubing connector and 2 port needles	1 unit, delivered sterile
ATS5061		Tubing Connector Connection piece made from titanium for catheter-to-catheter connection, e.g. for port-revision with ATOMS implant. Includes catheter spare piece of 200 mm	2 units in 1 box, delivered sterile
ATS5031		Tubing Plug Made from titanium For temporary closure of ATOMS catheter	2 units in 1 box, delivered sterile
SFN 0930 G		Port Needle Port needle for adjustment for the ATOMS Scrotal Port Needle Length: 30 mm, 20 G x 30 mm	25 needles / box, delivered sterile

Urology Demo Products


Order Code		Product description	Included in delivery
ATS5041-DEMO		A.M.I. Demo ATOMS ATOMS implant	1 unit, delivered non-sterile
ATS-BONE		ATOMS Male Pelvis Demonstration Model Pelvic bone with ATOMS implant	1 unit, delivered non-sterile
ATS-TUNNEL		ATOMS Tunnelling Demo Pelvic bone model with soft rubber inner parts for demonstration of correct tunnelling Tunnellers not included	1 unit, delivered non-sterile

All demo products are strictly for demonstration purposes only and not intended for use in the operating theatre!

V20250306

A.M.I. Headquarters:

A.M.I. Agency for Medical Innovations GmbH
Im Letten 1
6800 Feldkirch
Austria
T +43 5522 90505 0
F +43 5522 90505 4006
e info@ami.at



Empowering Innovation.
Empowering Lives.

www.ami.at