

MACROLUX

Healing for life, and life for healing



Series	HC29	HO35	HC35	HO42
Insertion tube outer diameter	2.9 mm	3.5 mm	3.5 mm	4.2 mm
Instrument channel width	/	1.8 mm	/	2.0 mm
Compatibility	/	5 Fr	/	5 Fr
Infusion	Single infusion	Continuous infusion		
Pixel	160 k		920 k	
Field of view	120 °		140 °	
Direction of view	15 °			
Depth of field	3-100 mm			
Working length	220 mm			
Weight	50 g			

Reference:
*Angelis, C. D. 2003. "Office Hysteroscopy and Compliance: Mini-Hysteroscopy versus Traditional Hysteroscopy in a Randomized Trial." Human Reproduction 18(11): 2441-45. doi:10.1093/humrep/deg463.

Healing for life, and life for healing

- Minimizing Environmental Impact
Our disposable endoscopes are made with environmental care.
- Caring for the Environment, Caring for You
Let's invest in a sustainable future together.



MACROLUX
Healing for life, and life for healing

301, Building 3, NamTai Inno Park In Guang Ming Avenue,
Guangming, Shenzhen, P.R.China
Tel:0755-88656680 Website:www.microlite.cn

MACROLUX

SeleneView


Single-Use Digital Hysteroscope ▼



H series

HC29 | HC35 | HO35 | HO42

Office | Outpatient | Ambulatory Surgery Room ▼



“See and Treat”
with single use hysteroscope



Lower Cost

Shorten hospital stay and reduce hospitalization costs for patients.



High Efficiency

Same-day surgery and discharge




Less Pain


Reduce cervical trauma and preserve fertility

Gentle Care ▼

Comfort in every moment




Vaginoscopic approach offers three key advantages:
No need for cervical dilation, no anesthesia, and no speculum.




Superior Visibility

Clarity in a new way ▼



Chromvue



HD Pixel with 140° Wide FOV

Optimized Efficiency

Convenient in every step ▼



Ultra light 01

Only 50 g weight, comfortable for whole day operation



See and treat 02

Compatible with 5 Fr instrument



Take and use 03

No need for assembling



5s Fast uterine distension

Achieve fast clear vision

04